

700X & 750X & 880E API 文档

本文文件是 700X & 750X & 880E 的应用程序开发接口(Application Programming Interface)。

文文件修改记录

版本	日期	修改内容	修改人
1.0	2019/09/04	初版	Sandy Lin
1.1	2019/12/27	修改 Samling rate 及相关程序代码说明	Sam Hsu

目录

常數定義 錯誤! 尚未定義書籤。

藍牙特徵值	錯誤! 尚未定義書籤。
藍芽指令	錯誤! 尚未定義書籤。
取得 User 1 或 User2 筆數	錯誤! 尚未定義書籤。
取得 User 1 或 User2 資料	錯誤! 尚未定義書籤。
取得 ECG rawdata	錯誤! 尚未定義書籤。
取得裝置資訊.....	錯誤! 尚未定義書籤。
同步裝置時間.....	錯誤! 尚未定義書籤。
回到待機模式(進行 Realtime 測量).....	錯誤! 尚未定義書籤。
Realtime 測量	錯誤! 尚未定義書籤。
測量結果	錯誤! 尚未定義書籤。
刪除 User 資料	錯誤! 尚未定義書籤。

常数定义

名称	值	说明
BT_HEADER	1	BLE 接收测量结果
BT_STANDBY	2	待机模式
BT_MEASURE	3	BLE 接收测量数据
BT_DOWNLOAD	4	下载
BT_DOWNLOAD_WAIT	0	下载结束进行等待
BT_DOWNLOAD_HEADER	2	下载 Device 记录信息
BT_DOWNLOAD_RAWD	3	下载 ECG rawdata
BT_DOWNLOAD_U1_COUNT	5	读取 User1 笔数
BT_DOWNLOAD_U2_COUNT	6	读取 User2 笔数
BT_SETUP	5	设定 Device (时间)
BT_SETUP_700X	2	设定 Device (时间)
BT_ERASE_FLASH	6	删除记录
BT_ERASE_USER1	3	删除 User1 记录
BT_ERASE_USER2	4	删除 User2 记录
BT_CONFIG_INFO	7	读取 Device 信息
BT_CONFIG_INFO_DEVICE	1	读取 Device ID
BT_CONFIG_INFO_SETTING	2	读取 Device FW_version
Realtime(ChannelNo)		
CH_HR	70	ECG 测量读取心率
CH_MMHG	71	BP 测量读取 mmhg
CH_ECG	73	ECG rawdata(心电图)目前 Realtime 125 个资料为 1 秒, 收到值为 109 为 1mV
CH_AUTOSCALE	74	调整放大 (心电图) 1、2 为一倍;3 为两倍;4 为 3 倍 (此功能非必要)

蓝牙特征值

Service UUID - Characteristic UUID

A000-2a36 数据传输(读

-2a37 数据传输(写

蓝芽指令

取得 User 1 或 User2 笔数 (880E 只有 User1)

取得 User 1 或 User 2 笔数

主要 Function - **public static void** download_user(**int** user)

function 传入参数:

- **user**: 取得数据之用户

BLE 传送指令:

// user = 1:

- **data2[1]** = **BT_DOWNLOAD**;
- **data2[2]** = **BT_DOWNLOAD_U1_COUNT**;

// user = 2:

- **data2[1]** = **BT_DOWNLOAD**;
- **data2[2]** = **BT_DOWNLOAD_U2_COUNT**;

mBluetoothLeService.writeCharacteristic(characteristic2, data2);

BLE 传回:

于 DataList functionBLE(**final int**[] value)

```
{
    format = value[1];
    bpCmd = value[2];
    if (format == BT_DOWNLOAD) {
        int recordsCount1, recordsCount2;
        if (bpCmd == BT_DOWNLOAD_U1_COUNT) {
            recordsCount1 = value[3];    //User 1 笔数
            recordsCount2 = 0;
        }
        else if (bpCmd == BT_DOWNLOAD_U2_COUNT) {
            recordsCount1 = 0;
            recordsCount2 = value[4];    //User 2 笔数
        }
    }
}
```

取得 User 1 或 User2 资料(880E 只有 User1)

取得 User 1 或 User 2 资料

主要 Function – 取得笔数后直接传送此指令

BLE 传送指令:

//判断数据笔数大于 0, 则传送取得数据指令

```
data2[1] = BT_DOWNLOAD;
```

```
data2[2] = BT_DOWNLOAD_HEADER;
```

```
data2[3] = (byte) headerSeq; // headerSeq: 从 0 开始 为第一笔数据
```

```
mBluetoothLeService.writeCharacteristic(characteristic2, data2);
```

BLE 传回:

于 DataList functionBLE(final int[] value) {

```
    format = value[1];
```

```
    bpCmd = value[2];
```

```
    if (bpCmd == BT_DOWNLOAD_HEADER) {
```

```
        seq= value [4];
```

```
        year= value [5]; //年
```

```
        month= value [6]; //月
```

```
        day= value [7]; //日
```

```
        hour= value [8]; //时
```

```
        minute= value [9]; //分
```

```
        second= value [10]; //秒
```

```
        UserMode= value [11] //目前下载的 user
```

```
        HeartRate= value [13]//心率
```

```
        HighBloodPressure = value [14] + value [15] * 256; //SYS
```

```
        LowBloodPressure = value [16] + value [17] * 256; //DIA
```

```
        if (HighBloodPressure == 0 && LowBloodPressure == 0)
```

```
            AnalysisType = TYPE_ECG; //ECG 量测模式
```

```
        else
```

```
            AnalysisType = TYPE_BP; //BP 量测模式
```

```
        WHOIndicate= value [18];
```

```

        value [19]& 0x01 ==0x01 // 显示 OK

        value [19]& 0x02 ==0x02 // ECG noise , HR 显示'EE'

        value [19]& 0x04 ==0x04 // ECG rhythm

        value [19]& 0x08 ==0x08 //ECG wave

        value [19]& 0x10 ==0x10 // ECG pause

        value [19]& 0x20 ==0x20 // ECG 心率 fast

        value [19]& 0x40 ==0x40 // ECG 心率 slow

        value [19]& 0x80 ==0x80 // BP 有 AF

    }

}

```

取得 ECG rawdata

取得 ECG rawdata
<p>主要 Function – public static void download_file(RecordList header) {}</p>
<p>function 传入参数:</p> <ul style="list-style-type: none"> - header: 欲下载资料的 header <p>BLE 传送指令:</p> <pre> UserMode = (byte) header.UserMode; headerSeq = (byte) header.Seq; downloadBufSeq = 1; data2[1] = BT_DOWNLOAD; data2[2] = BT_DOWNLOAD_RAWD; data2[3] = UserMode; //要下载的 userMode, 0=user1, 1=user2 data2[4] = headerSeq; //要下载的 seq data2[5] = downloadBufSeq; //从 1 开始传值, 收到 256byte 传送下一个, 收满 12* 256 + 68 * 125 停止 mBluetoothLeService.writeCharacteristic(characteristic2, data2); </pre> <p>BLE 传回:</p> <pre> 于 DataList functionBLE(final int[] value) { format = value[1]; // BT_DOWNLOAD bpCmd = value[2]; // BT_DOWNLOAD_RAWD if (bpCmd == BT_DOWNLOAD_RAWD) { </pre>

```

int Cmd = value[3];

if (Cmd < 15) {

    for (int i = 0; i < 16; i++)

        FlashBuffer[Cmd * 16 + i] = (byte) value[i + 4];

}

else if (Cmd == 15) {

    for (int i = 0; i < 16; i++)

        FlashBuffer[Cmd * 16 + i] = (byte) value[i + 4];

    if (currentSize < 12 * 256 + 68 * 125) {

        for (int i = 0; i < 256; i++) {

            rawDataBuf[currentSize++] = FlashBuffer[i];

        }

    }

    if (currentSize >= 12*256+68*125) {

        data2[1] = BT_DOWNLOAD;

        data2[2] = BT_DOWNLOAD_WAIT;

        mBluetoothLeService.writeCharacteristic(characteristic2, data2);

        //下载 Rawdata 完成，转换成 ECG 数据

        for (int i = 0, j = 0; i < (12 * 256 + 68 * 125) / 2 ; i++, j += 2) {

            iFirstByte = (short) (0x00FF & ((short) rawDataBuf[j]));

            iSecondByte = (short) (0x00FF & ((short) rawDataBuf[j + 1]));

            rawData[i] = (short) (iFirstByte << 8 | iSecondByte);

        }

    }

    else {

        // 传送下一个

        data2 [1] = BT_DOWNLOAD;

        data2 [2] = BT_DOWNLOAD_RAWD;

        data2 [3] = (byte) UserMode;

        data2 [4] = (byte) Seq;

```



```

        data2 [5] = (byte) ++downloadBufSeq;

        mBluetoothLeService.writeCharacteristic(characteristic2, data);

    }

}

}

```

取得装置信息

取得装置信息
<p>主要 Function – public static void get_info() {}</p> <p>BLE 传送指令:</p> <p>//取得 DeviceID</p> <pre> data[1] = BT_CONFIG_INFO; data[2] = BT_CONFIG_INFO_DEVICE; mBluetoothLeService.writeCharacteristic(characteristic2, data); </pre> <p>//取得 FirmwareVersion</p> <pre> data[1] = BT_CONFIG_INFO; data[2] = BT_CONFIG_INFO_SETTING; mBluetoothLeService.writeCharacteristic(characteristic2, data); </pre> <p>BLE 传回:</p> <p>于 DataList functionBLE(final int[] value) {</p> <pre> format= value [1]; // BT_CONFIG_INFO; bpCmd= value [2] ; if (format == BT_CONFIG_INFO) { if (bpCmd == BT_CONFIG_INFO_DEVICE) { int Device ID = value[12] * 256 * 256 * 256 + value [13] * 256 * 256 + value [14] * 256 + value [15]; String DID= String.format("%08x", DeviceID) //继续取得 Firmware Version data[1] = BT_CONFIG_INFO; data[2] = BT_CONFIG_INFO_SETTING; </pre>

```

        mBluetoothLeService.writeCharacteristic(characteristic2, data);

    }

    else if (bpCmd == BT_CONFIG_INFO_SETTING) {

        FirmwareVersion = (value[3] * 256 + value[4]);

    }

}

```

同步装置时间

同步装置时间
<p>主要 Function – public static void time_setting() {}</p> <p>BLE 传送指令:</p> <pre> data[1] = <i>BT_SETUP</i>; data[2] = <i>BT_SETUP_700X</i>; data[3] = year; //年, 2019 年则设定为 19 data[4] = month; //月 data[5] = day; //日 data[6] = hour; //时 data[7] = minute; //分 data[8] = second; //秒 mBluetoothLeService.writeCharacteristic(characteristic2, data); </pre> <p>BLE 传回:</p> <p>于 DataList functionBLE(final int[] value) {</p> <pre> format= value[1]; // BT_SETUP; if(format = <i>BT_SETUP</i>) { /" 时间设定成功"/ } } </pre>

回到待机模式(进行 Realtime 测量)

回到待机模式(Realtime)
主要 Function – public static void standby() {}
BLE 传送指令: data[1] = BT_STANDBY mBluetoothLeService.writeCharacteristic(characteristic2, data);

Realtime 测量

Realtime 测量
BLE 传回: 于 DataList functionBLE(final int[] value) { format = value [1]; if (format == BT_MEASURE) { int j = 0; for (int i = 4; i < 20;) { ChannelNo = value[i++] & 0xff; ChannelMSB = value[i++] & 0xff; ChannelData = (value[i++] & 0xff) * 256; ChannelData = ChannelData + (value[i++] & 0xff); if (ChannelNo == CH_HR) { //HR HeartRate = ChannelMSB * 256 * 256 + ChannelData; } else if (ChannelNo == CH_MMHG) { //BP bpDiastolic = ChannelMSB * 256 * 256 + ChannelData; } else if (ChannelNo == CH_AUTOSCALE) { ecgSize = ChannelMSB * 256 * 256 + ChannelData; } else if (ChannelNo == CH_ECG) { //ecg data MDRawData = ChannelMSB * 256 * 256 + ChannelData; if (ecgCount < ecg_rawData.length) { //125*34, 共 34 秒 ecg_rawData [ecgCount] = (short) (MDRawData); ecgCount ++; if (ecgCount > 125 * 3) { //从第 3 秒开始画

```
        displayData[displayCount] = ecg_rawData[ecgCount - 1];
        displayCount++;
    }
}
}
}
```

测量结果

测量结果
<p>BLE 传回:</p> <pre>于 DataList functionBLE(final int[] value) { format= value [1]; if (format == <i>BT_HEADER</i>) { seq= value [4]; year= value [5]; //年 month= value [6]; //月 day= value [7]; //日 hour= value [8]; //时 minute= value [9]; //分 second= value [10]; //秒 UserMode= value [11]//目前的 user HeartRate= value [13]//心率 if (value [14] != 0 && value [16] != 0) AnalysisType = TYPE_BP; //血压量测模式 else AnalysisType = TYPE_ECG; //ECG 量测模式 HighBloodPressure = value [14] + value [15] * 256; //SYS HighBloodPressure = HighBloodPressure & 0x00FF //转正十六进制数字显示 LowBloodPressure = value [16] + value [17] * 256; //DIA LowBloodPressure = HighBloodPressure & 0x00FF //转正十六进制数字显示</pre>

```

        WHOIndicate= value [18];

        value [19]& 0x01 ==0x01 // 显示 OK

        value [19]& 0x02 ==0x02 // ECG noise , HR 显示'EE'

        value [19]& 0x04 ==0x04 // ECG rhythm

        value [19]& 0x08 ==0x08 //ECG wave

        value [19]& 0x10 ==0x10 // ECG pause

        value [19]& 0x20 ==0x20 // ECG 心率 fast

        value [19]& 0x40 ==0x40 // ECG 心率 slow

        value [19]& 0x80 ==0x80 // BP 有 AF

    }

}

```

删除 User 数据(880E 只有 User1)

删除 User 数据
<p>主要 Function – public static void delete_user(int user) {}</p>
<p>BLE 传送指令:</p> <pre> /*删除 User1 数据*/ strArray[1] = BT_ERASE_FLASH; strArray[2] = BT_ERASE_USER1; //删除 User 1 数据 /*删除 User2 数据*/ strArray[1] = BT_ERASE_FLASH; strArray[2] = BT_ERASE_USER2; //删除 User 2 数据 </pre> <p>BLE 传回:</p> <p>于 DataList functionBLE(final int[] value) {</p> <pre> format= value [1]; if (format == BT_ERASE_FLASH) { //删除成功 } } </pre>

